


How do you know if equipment is safe to use at Mon Health Medical Center?

■ Medical equipment

- Inspection tags are current.
- Equipment is cleaned before patient use.
- Equipment is working properly (red tags are available through the Facilities Department to label broken equipment). Contact BioMed if medical equipment is broken (304-598-1439).



■ Other equipment?

- For appliances, ensure it has been electrical checked by Facilities.
- IT equipment should be issued or checked by IT.
- For utility failures, contact Facilities.

What to do in the event of an Environment of Care Incident?

- There is a Hospital Systems Failure and Basic Staff Response Quick Reference Guide attached to the Utility Systems Management Plan in PolicyStat.

- Snapshot of Guide:

HOSPITAL FAILURE	WHAT HAPPENS	WHO TO CONTACT	RESPONSIBILITY OF STAFF
Cardiac Monitoring System	Monitor Fails	BioMed x1439 (304-598-1439) House Supervisor x6648 (304-285-6648)	Use portable monitors on critical patients. BioMed will contact vendor for immediate service.
Computer Network	Computer system unavailable.	Call Helpdesk x1327 (304-598-1327)	Implement downtime procedures and remain off computer until notified.
Drinking Water	Contamination of Water Supply	Facilities x6680 (304-285-6680) Security x6664 (304-285-6664) House Supervisor x6648 Administrator on Call Incident Command All Department Managers	Facilities shuts off water supply from the outside until declared safe for use. Alert staff about water contamination and Do Not Use water for drinking. Facilities will call municipal utility company. Do not DRINK water; however, water can be used for other reasons. See <i>Handling Contamination of the Water Supply</i> policy and attachment.
Water	Loss of Water	Facilities Management x6680	Facilities Management and Security will initiate fire watch. Utilize bottled water and alcohol hand sanitizer. Facilities will call vendor for immediate service.
Electrical Power Failure – Generator starts	Only lighting and equipment on emergency power will operate. Only	Facilities Management x6680	Verify life support equipment is functioning. Use flashlights. Complete surgery cases in progress.